

# PixelSensor™ from PIXELTEQ

開發攜帶簡易型 — 最佳幫手！

## 規 格



積分時間：1-1024ms

電流：20-5120nA

8通道：

425,455,485,515,555,615,660,695nm


優點：價格便宜



425 455 485 515 555 615 660 695 (nm)

### •磁珠抗體螢光法

### •PCR/ RTPCR





<http://pic.biodiscover.com/topic/BEAVER/method.html>

<http://www.glass.cn>

•Sybr green 1倍，0.5倍  
0.25倍，0.125倍吸收值

•奈米金0.38倍，0.36倍，0.32倍，0.28倍  
0.25倍，0.22倍吸收值測試



Labguide Co., Ltd. (Taiwan)  
E-mail : terry@seed.net.tw  
Tel. : +886-2-29406162  
Add. : 3F., No.37, Ln. 232, Huaxin St., Zhonghe Dist., New Taipei City 235, Taiwan (R.O.C.)  
<http://www.labguide.com.tw>

Labguide Co., Ltd. (Shenzhen)  
E-mail : admin@labguide.com.cn  
Tel. : +86-075527750787  
Add. : Address: C603, Baohe Building, No.98 of Baonong North Alley, Bao'an District, Shenzhen, China.  
<http://www.labguide.com.cn>