

PRODUCT PORTFOLIO

GBC's Product Lines

AAS



HPLC



ICP-OES



UV-Vis



ICP- α TOFMS



Rheometry



XRD



4 Lakewood Boulevard
Braeside Victoria 3195
Australia

Telephone: +61 3 9588 6666

Fax: +61 3 9588 6677

Email: gbc@gbcsci.com

Web: www.gbcsci.com





XplorAA

- Wavelength range of 175 – 900 nm
- Continuously adjustable slit width from 0.1 to 2.0 nm for both normal and reduced height
- Extended D₂ working range of 175 – 425 nm
- Asymmetric modulation – 40% noise reduction
- HYPER-PULSE background correction for accurate correction of transient background signals to 2.5 Abs
- Choice of Single beam or Double beam optics
- Choice of manual 2 lamp or motorised 6 lamp turret
- Choice of Manual, Automatic or Programmable Gasbox
- Optional Super Lamp power supply for improved sensitivity, detection limits and calibration linearity
- Optional coded hollow cathode lamp recognition
- Optional Air Purge Kit which allows clean air from air cylinder or an air compressor to prevent acidic fumes entering the instrument



SavantAA Σ

- World's best performing AAS
- Wavelength range of 175 – 900 nm
- Continuously adjustable slit width from 0.1 to 2.0 nm for both normal and reduced height
- Extended D₂ working range of 175 – 425 nm
- Asymmetric modulation – 40% noise reduction
- HYPER-PULSE background correction for accurate correction of transient background signals to 2.5 Abs
- Automatic Burner Rotation for measurement of over range samples without the need for dilution
- Automatic Work Head Adjustment in vertical and horizontal direction for optimum sensitivity for all elements in a multi-element analysis
- Super Lamp power supply for improved sensitivity, detection limits and calibration linearity
- Electronic Sample Viewing (ESV) – enables view of furnace injection, sample drying and ashing from computer monitor to facilitate method development
- Coded hollow cathode lamp recognition



SavantAA Z enduro

- Wavelength range of 175 – 900 nm
- Continuously adjustable slit width from 0.1 to 2.0 nm for both normal and reduced height
- Variable Magnetic Field Strength from 0.6 to 1.1 Tesla in 0.1 Tesla increments
- 1 ms delay between pulses and measurements ensures best accuracy
- Electronic Sample Viewing (ESV) – enables view of furnace injection, sample drying and ashing from computer monitor to facilitate method development
- Coded hollow cathode lamp recognition



Cintra UV-Vis

- World's best specifications
- Wavelength range of 190 – 1,200 nm
- Choice of fixed or variable slit width
- Complete range of accessories
- 7 x 7 carousel allows 7 samples and 7 reference samples OR can be programmed for 12 samples and 2 reference samples for increased batch size sample analysis



OptiMass 9600 TOF ICP-MS

- World's only SIMULTANEOUS Time of Flight ICP-MS
- Mass range 1 to 260 amu
- Very low ppt detection limits
- 5 times faster sample throughput compared to Quadrupole ICP-MS
- Fast sample screening for all masses
- Isotope ratio precision comparable with TIMS
- Patented Smartgate for removal of interfering masses
- Octopole Collision Cell with mass independence for interference management during sample analysis



Emma X-ray Diffractometer

- The enhanced Multi-materials Analyzer
- 3 kW high power X-ray Generator
- Optional High speed multi-strip Solid State Detector (SSD)
- SSD Scan speed of 120 deg/min
- Theta/2Theta and Theta/Theta configuration
- High temperature and environmental stages
- Large and bulky sample stages
- Unique harmonic gearbox with zero backlash



Quantima ICP-OES

- 160 – 800 nm wavelength range
- 1800 / 2400 l/mm holographic grating
- Inbuilt water cooling
- 40 MHz free running RF generator allows fastest response to plasma changes
- Direct aspiration of organic and oil samples without any need for oxygen accessories
- Plasma Camera for sample viewing
- Patented PIP – significantly lowers occurrence of torch meltdowns
- Fully thermostatted optics enables excellent long term stability



ISO 9001
Quality Accreditation

GBC Scientific Equipment

Manufacturer of world class scientific equipment and accessories; AAS, HPLC, ICP-OES, ICP-oTOFMS, Rheometry, UV-Vis and XRD - since 1978.

Visionary Technology

GBC Scientific Equipment will advance people's knowledge and their capacity to enhance the quality of life for all humankind.

